

Primeprint Solution Checklist

As of August 2022





Dear Consultant, Dear Service Technicians,

This checklist assists you in preparing a quotation and during the smooth installation of Primeprint and Primeprint PPU to achieve the best possible customer satisfaction.

The Primeprint Solution Checklist:

- points out relevant technical requirements for a comprehensive quotation preparation
- shows the necessary technical logistic framework for a smooth installation
- provides basic information on the most important elements and functions of Primeprint Solution

Please note that this document does not replace the official Dentsply Sirona instructions for use for Primeprint, Primeprint PPU, software, and accessories.

We wish you success!

Your Dentsply Sirona Team

Contents

Prior to sale	04 - 05
Prior to installation	06 - 11
Day of installation	12 - 13
Cleaning and maintenance	14 - 19
Accessories and consumables	18 - 21

Prior to sale

In order to prepare a quote for Primeprint Solution, the following important parameters should be discussed and agreed with the customer. Please follow the recommendations below in order to ensure a smooth integration of Primeprint, Primeprint PPU, and the required accessories.

		Yes	No	Recommendation
PC	inLab PC* of at least version 5.0 or inLab 4 PC with performance package available			Include a quote for the suitable inLab PC for REF No. see pricelist
	System monitor available with a minimum resolution of 1920 x 1080, at 70 Hz			Include a quote for inLab system monitor REF 60 42 548
	Internet connection available			
Workstation	General Requirement			
	Set-up area is well ventilated**			
	Set-up area requirements for Primepri	nt		
	Loading capacity: min. 50 kg (110 lb)			
	Clearance height: min. 67 cm (26.37 in)			
	Flat and horizontal set-up area: min. $53 \times 51,5$ cm (20.86×20.27 in), ($W \times D$)			
	Space between the device and the wall of at least 5 cm (1,97 in)			
	Guaranteed access to the mains connector (AC 100 V - 240 V, 50/60 Hz, 2.0 A - 0.85 A)			
	Guaranteed access to network			
	Set-up area requirements for Primepri	nt PPU		
	Loading capacity: min. 60 kg (132 lb)			
	Clearance height: min. 67 cm (26.37 in)			
	Flat and horizontal set-up area: min. 73×55 cm (28.74 \times 20.27 in), (W \times D)			
	Space between the device and the wall of at least 8 cm (3,15 in)			
	Guaranteed access to the mains connector (AC 100 V - 240 V, 50/60 Hz, 4.2 - 2.2 A)			
	Guaranteed access to network			

If relevant, include further components in your quotation or contact other service providers together with your customer.

		Yes	No	Recommendation
Workstation	Set-up area for accessories			
	Space for at least one Primeprint Box 110 x 230 x 210 mm (4.33 x 9.06 x 8.27 in), (W x D x H)			
	Space for at least one Material Unit 280 x 190 x 120 mm (11.02 x 7.48 x 4.72 in), (W x D x H)			
	Space for at least two Washing Containers* each 110 x 230 x 210 mm (4.33 x 9.06 x 8.27 in), (W x D x H)			
Nitrogen	Space for a nitrogen bottle			
	Nitrogen with a purity of 2.6 which equals 99,6 Vol% available			If not: Contact a specialist seller and request a quote
	Additional accessories for the nitrogen supply:			If not: Contact a specialist seller and
	• Pressure regulator and manometer to regulate pressure to 4 to 8 bar			request a quote for a hose or an adapter that suits the current set up
	Primeprint PPU comes with a 5 m (196,85 in) hose (Ø 6 mm / 0,24 in), Push-in connector for 6 mm hoses with a G 1/4" inner thread, 5/16" x NPT ½" Push-in connector with an external thread, Push-In Reducer 5/16" x 6 mm and sealing cord. Depending on the individual set up additional accessories are needed.			
Cleaning agent	Isopropanol 99% available			If not: Contact a specialist seller and request a quote

^{*} inLab CAM Software is recommended to be installed on inLab PC.

^{**} For Primeprint Solution a dedicated fume cupboard is not required. A well ventilated room is generally recommended.

^{*} Two Washing Containers per material

Prior to installation

We recommend checking the conditions onsite prior to installation. This will ensure that everything runs smoothly on the day of installation. The most important points to consider are listed in the checklist below.

A. Contact list

r
ı



Customer				
Name:		_		
Street:				
Zipcode/town:				
Telephone:				
E-mail:				
B. Date of planned installation				
Date:	Time:	_		
If necessary post	onement to:			
Date:	Time:			



C. Structural requirements

Transport	Yes	No
Vehicle booked		
Transport for delivery to the installation site has been clarified. Packaging dimensions/weight: Primeprint: 830 x 620 x 1080 mm (32.68 x 24.41 x 42.52 in) (W x D x H) at 62 kg (136 lb) Primeprint PPU: 830 x 620 x 1080 mm (32.68 x 24.41 x 42.52 in) (W x D x H) at 72 kg (158 lb)		
Elevator available		
The necessary transport staff has been organized (2 persons)		

Set-up area	Yes	No
The minimum requirements for the Primeprint set-up area are fulfilled (see "prior to the sale")		
The minimum requirements for the Primeprint PPU set-up area are fulfilled (see "prior to the sale")		
The minimum space requirements for Primeprint Box, Material Unit and Washing Containers are fulfilled		

- A	

Connections	Yes	No
Pressure regulator available to reduce output pressure of the nitrogen to 4 to 8 bar		
Manometer available to check the supplied pressure of the nitrogen		
Hose with a diameter of 6mm or 6mm adapter is available		
Power sockets for PC, Primeprint and Primeprint PPU are available and freely accessible		
Network sockets for PC, Primeprint and Primeprint PPU are available and freely accessible		
Comments/measures:		



D. IT infrastructure

Hardware	Yes	No
PC with minimum hardware requirements available (see "prior to sale" and manufacturer information)		
System monitor available (see "prior to sale")		
Network	Yes	No
IT administrator is available on the day of installation (recommendation: Contact the IT administrator during the installation)		
Internet access for activating software licenses available		
Designated CAM PC is already integrated in the network		
Network with DHCP server available		
Automatic or static IP assignment in the network	Auto- matic	Static
TCP/IP address range		
Subnet mask		
Network connections with RJ45 connectors available close to each machine (min. 100 Mbit/sec; optimal 1 GBit/sec)		
Central storage location (shared network drive) available in network for restoration data		
Check availability of a PC or server that runs 24/7 (recommendation: Use an inLab PC or a server)		
Software	Yes	No
Ensure that all software requirements are met at the customer site. Observe the current manufacturer's specifications.		
Software type and version:		
Comments/measures:		

E. Equipment recommended on day of installation

Equipment	Yes	No
Trolley available		
Drill is available to drill holes to relocate wiring connections through worktops		
Test files (*.cam or *.stl) available for test printing and post processing		
Printing material is ordered and available on the day of installation		
Digital sliding caliper for measuring first print job		



Day of installation

On the day of installation, we recommend following the individual steps described in the Primeprint and Primeprint PPU user manual. In addition, go through the following installation protocol together with the customer. Ensure that printing material is available on site on the day of installation.

General	Yes	No
The user manual was handed over and the download of the instructions for use/software updates from the Dentsply Sirona homepage was mentioned		
inLab CAM Service was installed correctly and is running		
inLab CAM Software was installed correctly and licenses are activated		
Primeprint was installed as per the instructions for use		
Primeprint PPU was installed as per the instructions for use including securing and connecting the nitrogen bottle following the supplier information and loosening the transport lock		
The manometer shows a nitrogen pressure of 4 to 8 bar		
Primeprint can be started without any errors		
Primeprint PPU can be started without any errors		
The file exchange between the CAD Software and the CAM Software works flawlessly in the given folder/network structure		
An application was printed and post processed without errors (test cross, measurable)		

Briefing on inLab CAM key software functions	Yes	No
Creation of new production orders		
Import of restoration elements in *.CAM or *.STL format		
Selection and placing of elements on the building platform		
Positioning and adjustment of support structures		
Setting detail level		
Check CAM remote workflow		

Maintenance and machine parts	Yes	No
Explanation of the different materials and the color coding system		
Explanation of cartridges, Material Unit, Primeprint Box, and Washing Containers and briefing on their usage		
Briefing on changing the filters		
Briefing on cleaning agent handling		
Briefing on nitrogen handling		
Explanation of key error messages		
DS Core Care service & support*	Yes	No
Equipments registration on Siroforce (please be aware the service package will be available only for registered equipments)		
Explanation of service package		
DS Core Care flyer customer handover		
Conclusion	Yes	No
Contact information for ordering consumables and servicing was provided		



^{*} Service is subject to country availability.

Please contact your local Dentsply Sirona representative.

Cleaning and maintenance

Recommended maintenance intervals for Primeprint



Before every printing process

- Ensure material unit is locked
- Ensure cartridge is inserted and opened

Weekly

Check the glass pane of the light engine for resin residuals or other pollution and clean it if necessary. If the material unit is exchanged several times a week it is recommended to check it each time.

Monthly

- Exchange the inserted sponge in Primeprint Box
- Clean Primeprint Box with Isopropanol
- Clean the edge of the building platform with Isopropanol



Quarterly

 Check the collecting container for resin and if necessary clean it

Yearly

Annual maintenance carried out by a service technician

On request of the system

- Exchange the activated carbon filter
- Exchange the cartridge
- Exchange the foil in the vat

Before changing the print material

 Thoroughly clean the vat if it has been used for a different material before

Calibration

• The unit is calibrated at the factory

Recommended maintenance intervals for Primeprint PPU



Before post processing

- Ensure washing containers are unlocked
- Ensure that the cleaning agent fitness shown on the display is sufficient
- Ensure nitrogen bottle is opened and connected to the system



Weekly

• Check the ease of movement of the float magnet



Monthly

Check the cover of the light chamber for resin residuals

Yearly

Annual maintenance carried out by a service technician

On request of the system

- Exchange the activated carbon filter
- Exchange the ozone filter
- Exchange cleaning agent
- Exchange the nitrogen bottle
- Clean the washing containers

When changing the print material

• Exchange washing containers

Calibration

• The unit is calibrated at the factory

Tips for cleaning and maintenance



Clean the glass pane with Isopropanol if it is contaminated with resin, with lens cleaner (Purosol) and a microfiber cloth.





The level of the cleaning agent in the washing containers must always be between min and max.



The washing container can be locked (1) or unlocked (2). It is recommended to lock the washing containers when they are not used to reduce the loss of cleaning agent by evaporation.





The Primeprint
Solution material
cartridge comes with
colored stickers that
can be applied to the
handles of the washing
containers and to the
left side of the material
unit vat. This provides
a visual assignment
of the vat and the
washing containers.

Observe the local regulations for the disposal of all operating materials generated during the use of the Primeprint 3D printer such as cartridges, gloves or cleaning utensils.

Primeprint Solution

Accessories and Consumables

Product		REF.
	Primeprint Activated Carbon Filter	67 47 823
	Primeprint PPU Activated Carbon Filter	67 45 603
1500 E20	Primeprint PPU Ozone Filter	67 52 138
	Primeprint Box	67 44 895
	Primeprint Washing Container	67 45 546
	Primeprint Material Unit	67 44 903
	Transport Container Insert, 5 pieces	67 54 241
	Foil with Clamp Frame, Repl.	67 44 911

Primeprint Solution

Accessories and Consumables

Product		REF.
	Platform holder	67 56 402
	Side Cutter	67 58 663
Dentuphy •	Spatula DS	67 58 655
	Individual Curing Carrier	67 45 074
	Collecting Container, Repl.	67 57 434
	Microfiber Cloth	68 02 917
	Purosol	67 82 481
S Canada.	Playing cards	67 53 631

Primeprint Solution Material

Product		REF.	Description	Medical device Class MDR	Medical device Class FDA
	Primeprint Tray	67 40 216	Resin for individual impression trays	I	I
	Primeprint Guide	67 40 232	Resin for surgical guides	lla	1
	Primeprint Splint	67 40 240	Resin for hard splints	lla	1
	Primeprint Model	67 40 257	Resin for working models, situation models and control model	TEC resin	TEC resin
	Primeprint Model T	67 40 224	Resin for thermoforming models	TEC resin	TEC resin
	Primeprint Cast	67 40 265	Resin for dental casting objects for precision casting	TEC resin	TEC resin
	Primeprint Temp A1	67 40 281	Resin for temporary anterior and posterior tooth restorations A1	lla	Ш
	Primeprint Temp A2	67 40 299	Resin for temporary anterior and posterior tooth restorations A2	lla	II
	Primeprint Temp A3	67 40 307	Resin for temporary anterior and posterior tooth restorations A3	lla	Ш

Notes		
	-	
	-	
	-	
	-	
	-	
	-	
	-	
	-	
	-	
	-	
	-	

Dentsply Sirona

SIRONA Dental Systems GmbH Fabrikstraße 31, 64625 Bensheim, Germany dentsplysirona.com

